

Configuration Sheet – Shorts

Reference / Ordernumber: _____

(No Name!)

Step 1: What is ordered?

evomove orthokit (incl. holder)

right leg

female

Fabric color

evomove solokit (incl. pocket)

left leg

male


grey black

Pair of new shorts (excl. control unit)

Reference (Ordernumber, patient ID,..): _____

Step 2: Measurement (mandatory!)

Circumferences



U1

Waist
(midway between the lower costal arch and the iliac crest, level with the navel)

cm

U2

Buttocks
(strongest part buttocks/hips)

cm

Important:

Please measure both legs.

left

right

U3

Thigh
(1cm below crotch)

cm

cm


U4

Knee circumference
(back of the knee)

cm

cm

Lengths



Body dimensions are measured, e.g. using a **rubber band to mark** the positions **U1, U2** and **U4**.

L1

Knee length
(Waist to knee circumference U4)

cm

L2

Buttocks
(Waist – strongest point on the buttocks)

cm

L3

Seat length
(Waist to seat)

cm

Body height


cm

Important: The seat should be straight and without padding.

Step 3: Detection of the electrode position (please check)

○

1. Option: 3D-Scan



Create a **3D scan** using a scanner.

We accept **.obj** and **.stl** files.

Important: When scanning, please make sure that the buttocks and both legs can be seen standing from knee to waist.

○

2. Option: Photos




Photo with corresponding electrodes.

Important: Both legs must be visible from waist to knee when photographed.

Do not forget: Photograph laterally from both sides.

Step 4: Details of the shorts

1

Cord in waistband

- No
 Yes



2

Handypocket

- none
 2a right
 2b left

3

Pocket for evomove® control unit

- 3a right | medial | lateral
 3b left | medial | lateral



4

Zipper
(Desired position may be adjusted depending on the electrode position)

- No
 Yes, in the waistband
- 4a right | short | long
 4b center | short | long
 4c left | short | long

- Yes, on the leg
- 4d right | short | long
 4e left | short | long

Cable outlet
(standard in the pocket of the evomove® control unit)

- individual
(please highlight in the picture)

Wearing height

- Waist shorts | Hipster shorts